

COMPRESSOR
Aeras
 Ramvac

PREINSTALLATION GUIDE


MODELS:
 C3-2, C4-1, C4-2,
 C6-2, C7-2 & C11-2



Doctor: Address: Phone No.: Dealer: Address:	<hr/> <hr/> <hr/> <hr/> <hr/>
--	-------------------------------

Check <input checked="" type="checkbox"/> Compressor Model Being Installed:	<input type="checkbox"/> Aeras 3 (C3-2) <input type="checkbox"/> Aeras 4 (C4-1) <input type="checkbox"/> Aeras 4 (C4-2) <input type="checkbox"/> Aeras 6 (C6-2) <input type="checkbox"/> Aeras 7 (C7-2) <input type="checkbox"/> Aeras 11 (C11-2)
	<p>⚠ ATTENTION: It is recommended that this equipment be installed by a licensed electrical/plumbing contractor only. All installations must conform to local and national codes.</p>

NOTICE

 This information is not intended to replace the information found in the **User Manual**. Refer to the **User Manual** for full installation, safety and precautions.

PREINSTALLATION GUIDE




ELECTRICAL

HIGH VOLTAGE

- Install a dedicated electrical circuit of sufficient capacity (see [SPECIFICATIONS & SITE REQUIREMENTS](#) table).
- Any means provided to isolate this device from the supply mains shall isolate all poles simultaneously. The disconnection means and overcurrent protection are to be provided by the installer. This device should always be attached to a dedicated circuit with appropriate wiring and circuit protection.

LOW-VOLTAGE REMOTE SWITCHING (RECOMMENDED)



- Maximum wire length for low-voltage, 18-gauge wire is 500 feet.



AERAS CONTROL PANEL



LOW-VOLTAGE REMOTE SWITCHING (ALTERNATIVE)




Illuminated Remote Switch

SPECIFICATIONS & SITE REQUIREMENTS

	AERAS 3 (C3-2)	AERAS 4 (C4-1)	AERAS 4 (C4-2)	AERAS 6 (C6-2)	AERAS 7 (C7-2)	AERAS 11 (C11-2)
Max Users	3	4	4	6	7	11
Voltage (AC)	230 (208 - 253)	115 (103 - 126)	230 (208 - 253)	230 (208 - 253)	230 (208 - 253)	230 (208 - 253)
Power Consumption (Amps)	10.5	23	11.5	16.3	21	32
Breaker Size (Amps)	20	30	20	20	30	40
Fusetron Size	FRN 15	FRN 25	FRN 15	FRN 20	FRN 25	FRN 40
Number of Motor Heads	1	2	2	2	2	3
HP Per Motor Head	2	1.25	1.25	1.25 (Left) 2 (Right)	2	2
Pressure	85 PSI/586 kPa to 115 PSI/793 kPa with ± 2% Total Error Band					
High-Pressure Mode*	105 PSI/724 kPa to 131 PSI/903 kPa					
Weight (lbs)	160	190	190	250	260	385
Dimensions W x D x H (in)	29 x 21 x 29	29 x 21 x 29	29 x 21 x 29	35 x 25 x 30	35 x 25 x 30	41 x 29 x 32
Dimensions (with Sound Cover)	33.3 x 29 x 34	33.3 x 29 x 34	33.3 x 29 x 34	35 x 29 x 34	35 x 29 x 34	41.1 x 33 x 36
Peak Airflow Output @80 psi (CFM)	6.0	6.6	6.6	9.3	12.0	18.1
Pump Up Cycle Time, 85 to 115 psi (Sec)	45	45	45	45	40	40
Tank Size (Gal)	12	12	12	20	20	30
Filtration (microns)	0.3 Coalescing; 0.01 Final Particulate					
Dryer Duty Cycle	100%					
Dryer Dew Point (PDP)	-40°F / -40°C					
dBA Levels (with/without Sound Cover)	65/73	70/78	70/78	72/80	76/82	78/84

NOTES: *High-pressure mode should only be used with a compatible milling machine, not for use in general dentistry.

ETL CLASSIFIED
 Conforms to:
 ANSI/AAMI STD ES60601
 Certified to:
 CAN/CSA STD C22.2 NO. 60601-1
 Intertek
 4007136

© 2020 DentalEZ Alabama, Inc. DentalEZ is a registered trademark and Aeras is a trademark of DentalEZ, Inc. Ramvac is a registered trademark of Ramvac Dental Products, Inc.

Ramvac
 DENTALEZ

DentalEZ
 2500 Highway 31 South
 Bay Minette, AL 36507

Phone: 866-DTE-INFO
 Fax: 251-937-0461
 dentalez.com

PN: 7565-058B
 December, 2020