

# Forest

## DENTALEZ

### ALERT SYMBOL KEY

The safety notices associated with the following should be given special attention when they appear in maintenance, operating and emergency procedures in the guide.

**⚠ WARNING** indicates that the personal safety of the patient, end user or technician could be compromised by disregarding the WARNING. Not following instructions may result in an injury.

**⚠ CAUTION** indicates that a particular procedure or precaution must be followed to avoid system irregularity or possible damage to the product.

**! NOTE** indicates special information to improve the ease of maintaining the product, or to clarify important information.

**⚠ WARNING:** We recommend consulting a licensed professional, I.E. Structural or Architectural Engineer when attaching a delivery to the building structure. If installing the delivery on a cabinet, shelf or other furniture alternative, take extra precaution to assure the apparatus is structurally sound and secured to the building structure to support the delivery weight preventing detachment or tipping resulting in injury. Forest™ is not liable for injury resulting from improper installation.

**⚠ CAUTION:** Field modifications that alter the electrical and/or mechanical safety of Forest products conflict with agency construction file requirements and are not sanctioned by Forest.

**⚠ CAUTION:** Recommended and provided mounting hardware may not be suitable for your specific installation. Please use caution when selecting hardware to avoid possible damage. Forest is not liable for damage due to hardware choice.

**! NOTE:** Installation Verification: Refer to Operator's Guide and perform all operations. If operations performing as intended installation complete. If operations not performing as intended, review steps in installation instructions to confirm correct installation or call Technical Support at 800-423-3555.


Tools Required: Flathead and Phillips Screwdrivers.

### CABINET MOUNTING

STEP 1: Use product or dimensions below to install delivery (FIG. 1).

**! NOTE:** Depending on counter type, use appropriate screws to install unit (mounting hardware not included).

FIG. 1:



**! NOTE:** Not to scale.

# DELIVERY MOUNTING

STEP 1: Use product to mark delivery mounting screw placement (FIG. 1).

STEP 2: Secure unit under surface (FIG. 2).

STEP 3: Install clean water system bottle (FIG. 3).

FIG. 1:




FIG. 2:




FIG. 3:

