

Solara® QT High-Speed Handpieces

Lightweight and Powerful

- Titanium - 40% lighter than other Star® handpieces
- Neck designed specifically for smaller fingers.
- LubeFree ceramic bearings are more than 50% lighter than the bearings used by competitors.
- Engineered for more power when you need it most.
- Head size is smaller than competitor's head size for better visibility.

Ordering Information

Solara® QT Handpiece, Fiber Optic, LubeFree – #265718

Solara® QT Handpiece, Fiber Optic, Lubricated – #265791



StarDental®
Made in the U.S.A.



Comfort, Control and Value

Our StarSimplicity handpieces provide exceptional value and performance. The durable construction and thoughtful design allow precise control and greater efficiency in the dental operator.

- Durable stainless steel construction
- Lubricated ceramic bearings turbine for smooth, high-torque operation
- Improved unique water spray flushes site and keeps tooth and bur cool



StarSimplicity SW – #265633

Non-fiber-optic high-speed handpiece with lubricated ceramic bearing autochuck turbine. Uses standard 4-line StarDental swivel. 1-year limited warranty.



StarSimplicity FX – #265634

High-speed 4-line fixed back-end handpiece with lubricated ceramic bearing autochuck turbine. 1-year limited warranty.



Personalize your handpiece with your group name or logo!



For more information, please visit www.DentalEZ.com