

430 Series High-Speed Handpieces

The Power of Performance

The 430 series of handpieces from StarDental delivers the right combination of power and choice to drive your productivity. The 430 Torque has the highest power ratings tested at 27 watts*.

430 Series

The 430 Series handpieces are ergonomically designed with a stylish stainless steel or satin finish.

- Available in LubeFree, lubricated, fiber optic, and non-fiber optic
- **Small head design** — Size and shape promote superior operatory field visibility
- **Double bucket rotor design** — Delivers powerful performance
- **Vortex washer** — Unique washer designed to minimize debris “suck back” into housing, preventing turbine damage
- **Bifurcated glass** — Dual beam glass design prevents shadowing in front of bur (*F.O. models only*)
- **Co-axial water spray** — Optimizes spray, cooling the bur and flushing the operatory site
- **100% stainless steel body** — Rugged durable construction



Choose the Right 430 for Your Practice



430 Torque Satin Finish

High-speed handpiece that provides 27 watts of high torque power in a sleek, ergonomic design.

- 430SWL Torque LubeFree
- 430SWL Torque Lubricated
- 430SW Torque LubeFree
- 430SW Torque Lubricated



430SWL® Stainless or Satin Finish

Fiber optic high-speed handpiece features bifurcated glass to prevent shadowing from the bur.

- 430SWL LubeFree
- 430SWL Lubricated
- 430SWL Satin LubeFree



430SW™ Stainless or Satin Finish

Non-fiber optic high-speed handpiece with a rugged, 100% stainless steel body.

- 430SW LubeFree
- 430SW Lubricated
- 430SW Satin LubeFree
- 430SW Satin Lubricated



430 Flex and 430 Torque Flex Satin Finish

High-speed handpiece connects to a multiflex coupler and provides up to 31 watts of power.

- 430SWL Torque Flex LubeFree
- 430SWL Torque Flex Lubricated
- 430SWL Flex LubeFree
- 430SWL Flex Lubricated



430SWL 45 Stainless Finish

High-speed handpiece with a 45 degree head to allow better access to the posterior.

- 430SWL 45 LubeFree
- 430SW 45 LubeFree



430K Stainless Finish

High-speed handpiece features a fixed backend that does not require a swivel.

- 430K LubeFree
- 430K Lubricated

Personalize your handpiece with your group name or logo!

Ask your StarDental Customer Service Representative for more information.

LubeFree products include a 1-year limited warranty. Lubricated products include a 2-year limited warranty. Bifurcated solid cellular glass rod fiber optics include a 5-year limited warranty.



Please visit www.DentalEZ.com/StarDental for more information.