

## Titan® Low-Speed Handpieces

### Elite Comfort and Performance

Titan low-speed handpieces have a balanced, ergonomic design which helps to reduce operator fatigue. Titan motors, attachments and lubefree angles don't require lubrication and are adaptable to all your low-speed needs. All lubefree angles include a 1 year warranty.

### Titan Series

#### Titan 3 and Titan-T™

- Speed range 100-5,000 rpm (5k motor)
- Speed range 100 - 20,000 rpm (20k motor)
- Titan 3 has a timeless stainless steel finish
- Titan-T is made from titanium and offers a sleek, modern finish
- Exceptional balance
- 360° quick connect swivel

#### Titan E

- Lightweight aluminum exterior
- Reliable stainless steel E-coupling for attachment compatibility
- Available in fixed back end and swivel versions
- Fits all ISO standard e-type contras and attachments



Procedure	LubeFree System Configuration Options				Angles	Operating RPM Range		
<b>Prophy</b>	5,000 RPM Motor		Motor-to-Angle Adaptor	Five Star Prophy		1,000 - 3,000		
	20,000 RPM Motor	Torque Multiplier	Motor-to-Angle Adaptor	Five Star Prophy		1,000 - 3,000		
<b>Grinding &amp; Trimming</b> Denture, crown & bridgework, adjustment, disking, fabrication and trimming of temporaries	20,000 RPM Motor		Straight Nose Cone			10,000 - 20,000		
<b>Post &amp; Pin Drilling</b>	5,000 RPM Motor		Motor-to-Angle Adaptor	BB Manual Latch	BB Auto Latch	1,000 - 5,000		
	20,000 RPM Motor	Torque Multiplier	Motor-to-Angle Adaptor	BB Manual Latch	BB Auto Latch	1,000 - 5,000		
<b>Pin Setting</b>	5,000 RPM Motor		Motor-to-Angle Adaptor	BB Manual Latch	BB Auto Latch	1,000 - 5,000		
	20,000 RPM Motor		Motor-to-Angle Adaptor	10:1 Gear Reduction Latch		500 - 2,000		
	20,000 RPM Motor	Torque Multiplier	Motor-to-Angle Adaptor	10:1 Gear Reduction Latch		500		
<b>Finishing</b> Cavity and crown preps	20,000 RPM Motor		Motor-to-Angle Adaptor	BB Manual Chuck Friction Grip	BB Autochuck Friction Grip	BB Auto Latch	BB Manual Latch	10,000 - 20,000
<b>Polishing</b> Amalgam, composites, gold margins	20,000 RPM Motor		Motor-to-Angle Adaptor	BB Manual Latch	BB Auto Latch			10,000 - 20,000
<b>Endodontic</b> Canal enlargements (using Gates-Glidden drills) and sealing (using Lentulo), debridement (using file reamers) and reamers	5,000 RPM Motor		Motor-to-Angle Adaptor	Endodontic Auto Latch	Endodontic Manual Latch			350 - 5,000
	20,000 RPM Motor	Torque Multiplier	Motor-to-Angle Adaptor	10:1 Gear Reduction Latch				500 - 2,000
<b>Endodontic</b> Condensing gutta percha	20,000 RPM Motor		Motor-to-Angle Adaptor	BB Manual Latch	BB Auto Latch			10,000 - 20,000
<b>Endodontic</b> Root Canal Filing	5,000 RPM Motor		16:1 Contra Angle Attachment	BB Manual Latch	BB Auto Latch			300
	20,000 RPM Motor	Torque Multiplier	16:1 Contra Angle Attachment	BB Manual Latch	BB Auto Latch			300
<b>Caries Removal</b>	5,000 RPM Motor		Motor-to-Angle Adaptor	BB Manual Chuck Friction Grip	BB Autochuck Friction Grip	BB Auto Latch	BB Manual Latch	2,000 - 5,000
	20,000 RPM Motor	Torque Multiplier	Motor-to-Angle Adaptor	BB Manual Chuck Friction Grip	BB Autochuck Friction Grip	BB Auto Latch	BB Manual Latch	2,000 - 5,000



Please visit [www.DentalEZ.com](http://www.DentalEZ.com) for more information.