


ADAPTER INSTALLATION SHEET

(4 PIN SWIVEL TO 6 PIN TUBING)

VOLTAGE REQUIREMENT = 3.3VOLTS MAXIMUM



1. ALIGN AND THREAD TUBING
ASSEMBLY ONTO ADAPTER

2. ALIGN AND THREAD ADAPTER
ONTO TYPE-B SWIVEL

StarDental® Instrument Solutions

REPLACEMENT PART:
LAMP MODULE, REPLACEMENT
STARDENTAL PART NUMBER 265564