

COMPRESSORS

# Aeras

Ramvac

Aeras™ Compressors are smart enabled

Aeras Compressor Warranty 6 years/4,200 hours, extended to 8 years/5,600 when Smart Monitoring is enabled\*

## Electrical Requirements

Model Number	Description	Number of Users	Voltage	Dimensions W x D x H (in)	Break Size	Peak Airflow Output @ 80 PSI (CFM)
C3-2	Aeras 3 Compressor	3	230V	29" x 21" x 29"	20 amp	6.0
C4-1	Aeras 4 Compressor	4	115V	29" x 21" x 29"	30 amp	6.6
C4-2	Aeras 4 Compressor	4	230V	29" x 21" x 29"	20 amp	6.6
C6-2	Aeras 6 Compressor	6	230V	35" x 25" x 30"	20 amp	9.3
C7-2	Aeras 7 Compressor	7	230V	35" x 25" x 30"	30 amp	12.0
C11-2	Aeras 11 Compressor	11	230V	41" x 29" x 32"	40 amp	18.1

\*Terms apply

VACUUMS

# Aeras

Ramvac

Aeras™ Dry Vacuum Warranty 6 years, 8 years when Smart Monitoring is enabled\*, 10 years on rotary vane pump  
 Badger Warranty 2 years, 5 years on pump only.

## Specifications

Model Number	Description	Number of Users	Input Power	Breaker Size	Dimensions W x D x H (in)	Stacked Dim. W x D x H (in)
RV4-12T	Aeras 4	4	115v or 230v, 1 ph 60 Hz	15 amp	25" x 17" x 28"	29" x 19" x 48"
RV5-2T	Aeras 5	5	230v, 1 ph 60 Hz	20 amp	26" x 19" x 31"	33" x 21" x 59"
RV7-2T	Aeras 7	7	230v, 1 ph 60 Hz	20 amp	26" x 19" x 31"	33" x 21" x 59"
RV10-2T	Aeras 10	10	230v, 1 ph 60 Hz	20 amp	26" x 22" x 31"	33" x 22" x 59"
Badger LF**	Badger LF	4	208-230 VAC Single Phase	15 amp	25" x 18" x 24"	NA

\*Terms apply

\*\*Aeras monitoring not available with this product

©2021 DENTALEZ, Inc. Aeras is a trademark of DENTALEZ, Inc.  
 Ramvac is a registered trademark of Ramvac Dental Products, Inc.