

Star

400 FLEX

ADD STAR® POWER TO YOUR KAVO STYLE SWIVEL WITH SEAMLESS INTEGRATION

The latest addition to our handpiece family, Star 400 Flex with Aeras Intel, is compatible with the KaVo MULTIFlex™ coupler.

	Star 400 Flex	Competitor K*
Max cutting power supported by MULTIFLEX	YES	YES
Keep the power, lighten the load by 30%	YES	No
Exceptional ease of navigation via smaller head size	YES	No
RFID Technology embedded for powerful tracking and record-keeping**	YES	No
Made in the USA	YES	No

WHAT ELSE IS STANDARD WITH AN UPGRADE TO STAR?

Multi-port water spray

Cool-spectrum LED illumination

In-office turbine replacements

LubeFree Technology

Exceptional balance and in-hand stability

Star quality



*Competitive information based on Kavo model MASTERtorque LUXM8900L user manual.

** When used in conjunction with a compatible RFID reader system.

STAR HANDPIECES: SEAMLESS INTEGRATION WITH COMPETITOR SWIVELS



SPECIFICATIONS

STAR 400 FLEX

430SWL® TORQUE FLEX

430SWL FLEX

430SWL RA

Compatible with Kavo MULTIflex™ couplers	YES	YES	YES	No
Compatible with W&H couplers	No	No	No	YES
Weight	63 grams	60 grams	60 grams	51 grams
Head Size Diameter	11 mm	11 mm	11 mm	11 mm
Power	23 watts	20 watts	14 watts	16 watts
Operating Air Pressure	38-43 psi	40-42 psi	30-34 psi	44-46 psi
No Load Speed	436,000 rpm	340,000 rpm	430,000 rpm	430,000 rpm
Air Consumption	1.80 SCFM @ 43 psi	1.90 SCFM @ 40 psi	1.85 SCFM @ 34 psi	1.60 SCFM @ 45 psi
Turbine Warranty	LubeFree - 1 year Lubricated - 2 year	LubeFree - 1 year Lubricated - 2 year	LubeFree - 1 year Lubricated - 2 year	LubeFree - 1 year Lubricated - 2 year
Embedded RFID Tag	YES	No	No	No

ORDERING INFORMATION

Item	Description	Item	Description
266245	Star 400 Flex, Fiber Optic, LubeFree	266212	430SWL Torque Flex, Fiber Optic, Lubricated
266260	Star 400 Flex, Non-Fiber Optic, LubeFree	266153	430SWL Flex, Fiber Optic, LubeFree
266275	Star 400 Flex, Fiber Optic, Lubricated	266154	430SWL Flex, Fiber Optic, Lubricated
266276	Star 400 Flex, Non-Fiber Optic, Lubricated	266182	430SWL RA, Fiber Optic, LubeFree
266213	430SWL Torque Flex, Fiber Optic, LubeFree	266183	430SWL RA, Fiber Optic, Lubricated



Contact A Sales Rep

DENTALEZ
DENTALEZ.COM

© 2023 DENTALEZ, Inc. DENTALEZ, 430SWL, and Star are registered trademarks of DENTALEZ, Inc. MULTIflex is a trademark of the KAVO KERR GROUP. AERAS, Inc. is a separate corporation from DENTALEZ, Inc.