

INTERNATIONAL CATALOG 2022

# Affordable Excellence


Forest knows dentist fundamentally appreciate the value of many details done well and so do we. Our equipment is engineered with the highest quality materials so get the most from your investment, all for an affordable price.

# Designer Friendly®


Customize your equipment in Forest's exclusive Designer Friendly® finishes for a look that not only reflects your style, but catches the eye of your patients. Your equipment should be more than functional, it should elevate the value of your practice.

# Forest HANDS-FREE HVE

# Let's Dance.



Powder Coat Paint Color: Charcoal Ultraleather® Raven Wing & Sugar Powder Coat Paint Color: Charcoal Powder Coat Paint Color: Charcoal | Plastic Finish: Greywood Powder Coat Paint Color: Charcoal | Plastic Finish: Greywood



# TABLE OF CONTENTS





Chair Mounted Products & Configurator	08
Patient Chair	13
Delivery Control Heads	16
Mounting Options	18
Pivot	20
Fixed with Sidebox	22
Chair Mounted Lights + Posts	24
Chair Mounted Monitor Arms + Posts	25
Pre-Configured Chair Mounted Systems	26
Fixed Packages	27
Pivot Packages	31
7000 Series Rear Mount Side Delivery	36
Parts and Factory Options	38
	0.

# Navigating the Forest Catalog

The Forest catalog is organized into three sections:

Section

Chair Mounted Deliveries & Options

## Modular Chair Mounted Delivery Configurations

Use this section of the catalog to "custom" configure a chair mounted delivery package to your specifications. The product is designed to be "modular" starting with a chair simply add all desired products such as delivery control head, handpiece tubings, ancillary products, operatory lighting, monitor mounts, assistant's instrumentation or cuspidors.

- operatory chai
- chair mounted delivery control head
- pivot or fixed configurations, cuspidors & assistant's instrumentation
- chair mounted operatory lights
- chair mounted monitor mounts

Section 2

Pre-Configured Chair Mounted Packages

Forest has configured seven packages based on the most popular combinations of equipment with our pure white or Designer Friendly® color choice seletions

Section 3


7000 Series Rear Mount Side Delivery


# Step One

# Select A Patient Chair

Choose options for the selected patient chair

# 3900 Patient Chair







# Step Two 2

#### Select a Delivery Control Head

Choose desired handpiece tubing positioning & ancillary products such as electric micromotors, scalers, swivel trays or other options






# Step Three 3

# Select One of Two Mounting Options


Deliveries can be mounted on a chair with a left-to-right pivot mount or a fixed sidebox mount



# Step Four 4

#### Select Compatible Instrumentation & Cuspidors


Choose the assistant's instrumentation, cuspidors or related ancillary products




# Step Five

# Select Chair Mounted Lights or Monitor Display Mounts

Each light or monitor display requires a mounting arm; additionally the operatory light will require a power source if the 250 watt power supply was not selected in step two





# 3901INT



Shown with Memory Comfort Upholstery 3907-MC-WB

Model 3901INT (220V) \$9,550

Shown with Standard Asepsis Upholstery 3901INT (Plush upholstery is available at no additional charge.)

Standard Features	Factory Options	Part No	Price	Total
Ultraleather® upholstery with wide back	Hydraulic Chair (220V)	3901INT	\$9,550 _	
for elbow support	Hydraulic Chair (110V)	3900INT	\$9,550 _	
Double articulating pinch-style headrest	Designer Friendly® powder coat option	3901-220-CS	\$1,600 _	
Thin back with full lumbar support	Add chair touchpad	3914-102	\$650 _	
Right/left patient access	Less chair footswitch	3915-010-CS	-\$565 _	
Swing-down locking armrests	Quick disc. w/flow control 1/4" (for water)	0016-454-CS	\$190 _	
Hydraulic cantilever mechanism Emergency shut-off 30 degree base rotation with lock	Quick disc. w/flow control 3/8" (for air)  Chair adapter required for quick disconnects  Contoured utility center housing & cover	0016-456-CS	\$190 _	
Multifunction foot switch  w/4 programmable preset positions,	Standard with delivery head Contoured utility center housing & cover	3900-101	\$555 _	
accommodates 3 additional controls	Designer Friendly powder coat option	3900-101-C	\$675 _	
	Special color or upholstery	3907-999-SP	\$370 _	
Seat height: minimum 17" maximum 29"	Dual color upholstery	3907-999-DC	\$370 _	
220 Volt	Memory Comfort upholstery	3907-MC-WB	\$1,050 _	
	Soft-touch back massager	3900-100-INT	\$710 _	
	Chair Adapter comes standard with pivot and sidebox			
	For additional chair adapter selections see page 38			
			Total: _	

•	Designer Friendly® Surface Treatments								
	Powder coat paint for chair:								
	Standard colors	: O Forest Gre	ey O White						
	<b>OR</b> upgrade ch	air base/back p	oowder coat pa	int (add \$1,600):					
	O Arctic Silver	○ Bronze	$\bigcirc$ Charcoal						
	$\bigcirc$ Copper	$\bigcirc  Midnight$	$\bigcirc$ Pewter	○ Shadow					

۵	Ultraleather® Upholstery							
	Ultraleather® uph							
	O Autumn Leaf	○ Baltic		○ Ch	arcoal	○ Di	plomat Blue	
	○ Doe	O Dove Grey		○ Во	ysenberry	○ Pa	apyrus	
	○ Parrot	O Raven Win	ng	○ Sto	one		. John	
	Pearlized Colors:							
	○ Midnight	○ Pewter	$\circ$	)z	O Pharaoh	1	○ Tin Man	
	Special color (spe	ecify			_) 3907-999	9-SP	\$370	



# Standard IC Control Head

- Material: Aluminum hinged lid / 1/8" ABS plastic chassis with super tough DR acrylic protective coating
- Benefit: Hinge allows quick access to controls & ancillary products / lid & chassis are resistant to UV, impact / aluminum lid is strong & lightweight
- · Material: Electroless nickel plated brass air & water T's, routing valves, pneumatic control & distribution blocks
- Benefit: Lifetime warranty on all parts "under the hood" (on all Forest control heads (Max. 20 years) / trouble free handpiece operation
- · Material: Aluminum handpiece holder bar with glass filled nylon injection molded handpiece holders
- Benefit: Allows for up to four handpiece positions, a syringe & ancillary touch pad at your desired angle & spacing
- Material: Hinged aluminum control head lid
- Benefit: Available in nine two-coat powder coat paint colors
- · Material: Plastic chassis
- Benefit: Available in three carbon fiber look & two solid color finishes



# **Euro Control Head**

- Material: Whip arms rotate on nylon wheels with aluminum & brass components
- Benefit: Feather-touch whip movement / ultra light handpiece feel / long, trouble free life
- Material: Control head swivels at its center point
- Benefit: 240-degrees of control head rotation allows for beside or over-the-patient placement
- Material: Forward mounted aluminum whip arm system
- Benefit: Allows for greater reach of handpiece & syringe expanding control head placement options
- Material: Plastic control head lid
- Benefit: Available in three carbon fiber looks & two solid color finishes
- Feature: Autoclavable silicone handpiece pad available in black & white
- Benefit: Blends with Designer Friendly® choices

# **Delivery Control Heads**

#### 2001INT 2003INT



Standard IC Delivery
Model 2001INT \$5,275



Over the Patient Delivery

Model 2003INT \$6,740

#### Standard Features

Automatic Control For Up To Four Handpieces

For 2001INT & 2003INT

Asepsis 3-way air water syringe w/tubing, 5'

Contoured utility center, housing & cover

Pneumatic master control block w/diaphragm

Master on/off toggle

Individual handpiece pressure adjustment

Individual water coolant adjustment

Non-retraction water coolant with on/off

Air coolant adjustment

Handpiece exhaust filter

Single handle brake release (2001INT) dual with pivot mount

Dual handle brake release function (2003INT)

Wet/dry foot control

Air kit for utility center

Flexible control arm with pneumatic brake

Internal umbilical, 5' with 1' external from chair to utility center

\*Please note: tubing sold separately.

Shipping weight: 34 lb / 16 kg 42 lb / 19 kg

# Designer Friendly® Surface Treatments

Powder coat paint color:

Standard colors: O Forest Grey O White

**OR** upgrade powder coat paint (add 10%):

apgrado povidor code paine (add 1070).

○ Arctic Silver ○ Bronze ○ Charcoal

○ Copper ○ Midnight ○ Pewter ○ Shadow

#### Customize Delivery: Model Selected: Part No Price Total Step 1: Select delivery head model Standard IC: \$5,275 2001INT Step 2: Select surface treatment options Over-the-Patient: \$6,740 \_\_\_\_\_ 2003INT Step 3: Select handpiece positions

# Handpiece Positions:

Step 4: Select optional upgrades

Handpiece Positions	Description	1	2	3	4	5	Part No	Price	Total
SYRINGE	silicone sheath, grey	~					INCLUDED	N/C _	
4-HOLE (vinyl)	vinyl tubing, grey						0017-010-GRY	\$49 _	
■ 6-PIN (silicone)	silicone tubing, grey						0017-614-INT	\$250 _	
SCALER	Satelec SP Newtron Piezo scaler						SAT-S-INT	\$2,615 _	
LED SCALER	Satelec LED Piezo scaler						SAT-S-LED-INT	\$3,355 _	
s additional HP, wet	REQUIRED for a WET 4th HP position						1100-189-INT	\$540 _	
s ADDITIONAL HP, DRY	REQUIRED for a DRY 4th HP position						1100-190-INT	\$270 _	
NO TUBING	leave the block open						OPEN	_	
PLUG BLOCK	plug & close the block						PLUG	_	

REQUIRES HP illumination circuit board & <u>EITHER</u> 250 watt power supply or 25 watt power supply (board can accommodate (2) HP positions) See Page 40 when a fourth HP position is ordered, an ADDITIONAL DRY or WET HP position must be ordered


# Optional Upgrades:

Description	Part No	Price	Total
Upgrade Kit: wet/dry foot control with air chip	1106-002-INT-CS	\$100 _	
Unit mounted chair touchpad	1700-081-CS		
Additional umbilical (per foot)	1115-999-CS	\$45 _	
Water shut-off valve & regulator kit	1115-297-INT		
Quick disconnect w/flow control, 1/4" (for water)	0016-454-CS	\$190 _	
Quick disconnect w/flow control, 3/8" (for air)	0016-456-CS	\$190 _	
Chrome Accents	9000-998-CS	\$180 _	
HP illumination circuit board with time delay and adjustable voltage, accommodates 2 HP positions*			
*Requires power supply	0014-100	\$280 _	
Power supplies			
25 watt power supply (HP illumination)	0014-160	\$160 _	
48 watt power supply (motors)	0014-399	\$325 _	
250 watt (universal) power supply*	0014-387	\$680 _	
*Required when ordering two or more ancillary devices or when ordering chair mounted light			
Arm mounted swivel tray option			
Small Tray, side delivery (213mm x 318mm)	1110-468	\$440 _	
Small Tray, over the patient delivery (213mm x 318mm)	1110-787	\$415 _	
Large Tray, side delivery (289mm x 375mm)	1110-475	\$480 _	
Large Tray, over the patient delivery (289mm x 375mm)	1110-774	\$480 _	
Less contoured utility center (housing and cover only)	3900-101-CS	-\$250 _	
Add standard utility center (housing and cover only)	1115-020	\$190 _	
		Total: _	

# **Forest Mounting Options**


# Pivot Left/Right Mount

- Needle roller bearings & teflon friction adjustments 240 degrees of "feather-touch" rotation & maneuverability
- All aluminum arms, collars, cast mount adapter & water bottle mount creates a sturdy & strong foundation that will outlast the industry standard of plastic construction
- Chrome billet aluminum bearing covers & caps creates a rich & sturdy aesthetic



# Fixed With Sidebox Mount

- All aluminum housing with heavy duty mounting arm & feather touch bushings
- Large area to neatly store all ancillary products or house DURR CS1 amalgam separator
- High quality plastic trim, with billet aluminum caps, collars & plugs
- Regardless of the options you choose we provide a a customized feel for all option combinations
- All aluminum chassis, cast mount adapter & aluminum caps, plugs & collars available in nine two-coat powder coat paint colors
- Plastic trim & housing cover available in three carbon fiber looks & two solid color finishes



# 8601INT Pivot Configuration (Left/Right)



Model 8601INT \$2,055

#### Standard Features

Pivot chair adapter/ bearing kit Internal umbilical Pivot arm Single 1.5 liter self-contained clean water bottle Quick disconnect system for water bottle

**Delivery Head Options** 

The 8601INT pivot mount model must be ordered with one of the chair mounted delivery systems found in the delivery head configuration worksheets.

2001INT 2003INT

#### 

# Customize Delivery:

Step 1: Begin with selected mounting configuration

Step 2: Select surface treatment optionsStep 3: Select pivot configuration add-ons

Step 4: Select optional upgrades

Add-On	c.
- Auu-Oii	o.

Model

Price

Total



1095INT Assistant Instrumentation
with HVE/SE Syringe
Rear-mounted assistant's instrumentation
Quad-holder with HVE/SE and 3-way air/water syringe
Central vacuum accessory package with disposable solids collector

\$1,905 \_\_\_\_\_



1590INT Deluxe Cuspidor
Rear mounted cuspidor
Quad-holder with HVE/SE and 3-way air/water syringe
Central vacuum accessory package with disposable solids collector
One-piece ceramic bowl
Automatic cup fill and timed bowl rinse
Water shutoff valve and regulator kit

\$4,270 \_\_\_\_\_



1091INT-DURR Assistant's Instrumentation
Rear mounted assistant's instrumentation
with Durr electric-activated vacuum
Tri-holder with HVE/SE and 3-way air/water syringe
Central vacuum accessory package with disposable solids collector

\$2,405



1591INT-DURR Deluxe Cuspidor
Rear mounted cuspidor
with Durr electric-activated vacuum
Tri-holder with HVE/SE and 3-way air/water syringe
Central vacuum accessory package with disposable solids collector
One-piece ceramic bowl
Automatic cup fill and timed bowl rinse
Water shutoff valve and regulator kit

\$4,770 \_\_\_\_\_

# Optional Upgrades:

Description	Part No	Price	Total
Chair touchpad (1095 model only) Additional self-contained clean water system, single bottle 1.5 liter (1091 & 1095 models only)	1700-088-CS 1700-319	\$650 _ \$360 _ <b>Total</b> :	

# 8500INT Fixed With Sidebox Configuration



Model 8500INT \$2,915

#### Standard Features

Chair adapter

Sidebox housing

Internal umbilical

Control post

Bent rigid arm

Single 1.5 liter self-contained clean water bottle

Quick disconnect system for water bottle

**Delivery Head Options** 

The 8500INT fixed with sidebox model must be ordered with one of the chair mounted delivery systems found in the delivery head configuration worksheets.

2001INT 2003INT

#### 

## Customize Delivery:

Step 1: Begin with selected mounting configuration

Step 2: Select surface treatment options

Step 3: Select fixed with sidebox configuration add-ons

Step 4: Select optional upgrades

Α	ᅬ	٦		$\frown$	100	٠.
A	u	u	ᄆ	v.	401	S.

Model

Price

Total



1577INT Fixed Cuspidor Cuspidor (sidebox mounted) One-piece ceramic bowl Automatic cup fill and timed bowl rinse Water shutoff valve and regulator kit

\$2,375 \_\_\_\_\_



1576INT Deluxe Cuspidor
Cuspidor with assistant's instrumentation (sidebox mounted)
Quad-holder with HVE/SE and 3-way air/water syringe
Central vacuum accessory package with disposable solids collector
One-piece ceramic bowl
Automatic cup fill and timed bowl rinse
Water shutoff valve and regulator kit

\$3,195 \_\_\_\_\_



1092INT-DURR Assistant's Instrumentation
Assistant's instrumentation (sidebox mounted)
with Durr electric-activated vacuum
Tri-holder with HVE/SE and 3-way air/water syringe
Central vacuum accessory package with disposable solids collector

\$1,890 \_\_\_\_\_



1578INT-DURR Deluxe Cuspidor
Cuspidor (sidebox mounted)
with Durr electric-activated vacuum
Tri-holder with HVE/SE and 3-way air/water syringe
Central vacuum accessory package with disposable solids collector
One-piece ceramic bowl
Automatic cup fill and timed bowl rinse

\$3,695

## Optional Upgrades:

Description	Part No	Price	Total
Quick disconnect w/ flow control 1/4" (for water) Quick disconnect w/ flow control 3/8" (for air) Additional self-contained clean water system, single bottle 1.5 liter (1092INT model only) Durr amalgam separator cover (separator not included)	0016-454-CS 0016-456-CS 1700-319 1400-270-CS	\$190 _ \$190 _ \$360 _ \$500	
		Total:	

#### 9078INT 9083INT 9073INT





Standard Features



110V: available upon request LED Light Model 9078INT (220V) \$3,455



110V: available upon request LED Light Model 9073INT (220V) \$3,955

LED Light 9078INT
Automatic on/off activation with override
3-Axis light head movement
Easy release autoclavable handles
Light source 10 LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Power Supply, 100 Watt
Designer Friendly® powder coat paint light head

LED Light 9073INT
Automatic on/off activation w/override
3-Axis light head movement
Six-mode, polycarbonate button keyboard
Easy release autoclavable handles
Light source 8 1-watt LED lamps
Color temperature, 4,200 to 6,000 Kelvin
Designer Friendly powder coat paint light head

Shipping weight (est.): 35 lb / 16 kg

Factory Options	Part No	Price	Total
Light model choice:			
Light Post (For chair adaptors see page 38)			
Pivot (does not include cast chair adapter)	1113-654	\$760 _	
Pivot (includes cast chair adapter)	1113-654-CA	\$1,510 _	
35" Fixed w/clamp	1110-820	\$440 _	
44" Fixed two-piece light post	1110-616	\$365 _	
Short rigid arm 20" in lieu of standard 32"	1110-227-CS	\$N/C _	
Power Supplies			
The 100W power supply now comes standard wi	th INT lights		
Model 9083INT			



replacement bulb

L	Designer	⊢riendl∖	⁄™ Suri	tace I	reat	tments

\$52 \_\_\_\_\_

Total: \_\_\_\_\_

Powder coat paint color:
Standard colors: O Forest Grey O White

OR upgrade powder coat paint (add 10%):
Arctic Silver O Bronze Charcoal

○ Copper ○ Midnight ○ Pewter ○ Shadow

0014-300

## 1300RM-INT 1300TPM-INT





ring post mount Model 1300RM-INT \$1,260

top post mount Model 1300TPM-INT \$1,260

Standard Features	Factory Options	Part No	Price	Total
Monitor Mount Assembly, Ring Post Mount 1300RM-INT ICW Quicklink monitor mount Rex positioning handles 17.5" ring mount rigid arm	Monitor mount assembly model choice: Light Post (1300RM-INT) Pivot light posts & requirements 44" Fixed two-piece dual slot, sidebox mount w/ Forest light bushing	(call to order) 1110-617	\$365 .	
Monitor Mount Assembly, Top Post Mount 1300TPM-INT ICW Quicklink monitor mount Rex positioning handles 22" top post mount rigid arm	Monitor Post (1300TPM-INT) Pivot light posts & requirements 30" w/Dual slot, fixed sidebox mount, w/ Oilite bushing	(call to order) 1331-998	\$275 <sub>.</sub>	
Monitor Requirements: Weight: 8.25 to 13 lbs max, width of 19"			Total:	

All monitor arms must be ordered with a mounting post If ordering without a unit, a chair adapter is required with the purchase of the chair (monitors screens not included).

		Desigi	ner Friendl	y® Surfac	e Treatments
		Powder coat pa	int color:		
		Standard colors	s: O Forest Grey	/ O White	
Chinning weight:		<b>OR</b> upgrade po	wder coat paint	(add 10%):	
Shipping weight: 1300RM-INT, 1300TPM-INT	15 lb / 7 ka	O Arctic Silver	○ Bronze	○ Charcoal	
		<ul><li>Copper</li></ul>	<ul> <li>Midnight</li> </ul>	O Pewter	○ Shadow

Step One

Review The "As Shown Package"

Step Two

Review The Package Specifications, Part Numbers & Pricing

Step Three 3

Explore Additional Popular Equipment Choices Details; Can Be Found In "Factory Options & Parts" Located In Section 3



# Section 2 Pre-configured Chair Mounted Systems

\$30,647



Patient Chair		(page 13)	Operatory Lights + Monitor Mounts	(page 24)
Operatory chair	3901INT	\$9,550	*Halogen light 9083INT	\$2,790
Designer Friendly® Powder Coat	3901-220-CS	\$1,600	*44" Fixed two-piece light post	
Memory Comfort Upholstery	3907-MC-WB	\$1,050	with Dual Slot 1110-617	\$365
Dual Color Upholstery	3907-999-DC	\$370	*Ring Mount Monitor Assembly 1300RM-INT	\$1,260
			Subtotal	\$29,017
Delivery Head		(page 14)	*with 10% Designer Friendly upcharge	\$1,630
*Standard IC delivery	2001INT	\$5,275		
4-hole HP tubing, vinyl (quantity of 3)	0017-010-GRY	\$147	Total:	\$30,647
Configuration		(page 18)		
*Fixed with Sidebox package *Deluxe cuspidor	8500INT	\$2,915		
w/Durr assistant's instrumentation	1578INT-DURR	\$3,695		

# \$29,161



Patient Chairs Operatory chair	3901INT	(page 13) \$9,550
Delivery Systems Over the patient delivery 4-Hole vinyl HP tubing (quantity of 4) Small swivel tray (over the patient)	2003INT 0017-010-GRY 1110-787	(page 14) \$6,740 \$196 \$415
Configuration Sidebox package Deluxe cuspidor Durr assistant's instrumentation, Mounted from cast adapter	8500INT 1577INT 1092INT-DURR	(page 18) \$2,915 \$2,375 \$1,890

1110-617	\$365
1300RM-INT	\$1,260
9078INT	\$3,455
	(page 24)
	1300RM-INT

\$25,177



Patient Chair		(page 13)
Operatory chair	3901INT	\$9,550
Delivery Head		(page 14)
Standard IC delivery	2001INT	\$5,275
4-hole HP tubing, vinyl (quantity of 3)	0017-010-GRY	\$147
Small swivel tray (side delivery)	1110-468	\$440
Confinuation		( 10)
Configuration		(page 18)
Sidebox package	8500INT	\$2,915
Deluxe cuspidor		
w/Durr assistant's instrumentation	1578INT-DURR	\$3,695

Operatory Lights + Monitor Mounts		(page 24)
Halogen light	9083INT	\$2,790
44" Fixed two-piece light post	1110-616	\$365

Total: \$25,177

# Fixed With Sidebox

# \$26,666



Patient Chairs		(page 13)
Operatory chair	3901INT	\$9,550
Delivery Head		(page 14)
Over the patient delivery	2003INT	\$6,740
4-Hole vinyl HP tubing (quantity of 4)	0017-010-GRY	\$196
Large swivel tray (over the patient)	1110-787	\$415
Configuration		(page 18)
Sidebox package	8500INT	\$2,915
Deluxe cuspidor		
w/Durr assistant's instrumentation	1578INT-DURR	\$3,695

Operatory Lights + Monitor Mounts		(page 24)
Halogen light	9083INT	\$2,790
44" Fixed two-piece light post	1110-616	\$365

Total: \$26,666

\$23,422



Patient Chair Operatory chair	3901INT	(page 13) \$9,550
Delivery Head		(page 14)
Standard IC delivery	2001INT	\$5,275
Vinyl tubing 4-hole (quantity of 3)	0017-010-GRY	\$147
Small swivel tray (over the patient)	1110-468	\$440
Configuration		(page 18)
Pivot package	8601INT	\$2,055
Durr assistant's instrumentation	1091INT-DURR	\$2,405

Operatory Lights + Monitor Mounts		(page 24)
Halogen light	9083INT	\$2,790
Pivot light post	1113-654	\$760

Total: \$23,422

# Pivot (Left/Right)

\$26,452



Patient Chair Operatory chair	3901INT	(page 13) \$9,550
Delivery Head Standard IC delivery Vinyl HP tubing 4-hole (quantity of 3) Small swivel tray (side delivery)	2001INT 0017-010-GRY 1110-468	(page 14) \$5,275 \$147 \$440
Configuration Pivot package Pivot cuspidor	8601INT	(page 18) \$2,055
w/Durr assistant's instrumentation	1591INT-DURR	\$4,770

Operatory Lights + Monitor Mounts		(page 24)
LED light	9078INT	\$3,455
Pivot light post	1113-654	\$760

Total: \$26,452

\$27,341



Patient Chair		(page 13)
Operatory chair	3901INT	\$9,550
Delivery Head		(page 14)
Over the patient delivery	2003INT	\$6,740
Vinyl tubing 4-hole (quantity 4)	0017-010-GRY	\$196
Large swivel tray (over the patient)	1110-774	\$480
Configuration		(page 18)
Pivot package	8601INT	\$2,055
Pivot cuspidor		
w/Durr assistant's instrumentation	1591INT-DURR	\$4,770


Operatory Lights + Monitor Mounts		(page 24)
Halogen light	9083INT	\$2,790
Pivot light post	1113-654	\$760

Total: \$27,341

# Section 3 Customization

Step Two Step One Select Model & Mounting Position Select Surface Treatment Options Step Four Step Three Select Optional Upgrades Ancillary Products & Positioning

# **Cabinetry Mounting Configurations**



- 3/4" hardened steel 11"x 141/2" (pre-installed) mounting plate that bolts to the bottom of your counter top
- the 6"x 6" leveling plate is next, then a set of needle roller bearings with friction adjustment
- a super strong 21/2" steel tube scissor arm assembly
- · needle roller bearings are used where the two arm sections pivot

### 7000 Series: Rear Mount, Side Delivery



Model 7000INT-DURR \$8,770

Add Durr Vacuum Upgrade:



Optional storage organizer & swivel tray 1110-179-CS and 1110-909-CS





**Beneath Counter Mount** 7000-999-BC



**Under Cabinet Mount** 7000-999-UC



Wall Mount 7000-999-WM

#### Standard Features

#### **Duo Swing Delivery Head 7000 Series**

#### Duo control

automatic control w/lockout feature for three handpieces asepsis 3-way air water syringe w/tubing, 5' handpiece tubings factory installed at 4' (tubing is sold seperatly) pneumatic master control block w/diaphragm master on/off individual handpiece pressure adjustment

individual water coolant adjustment non-retraction water coolant

air coolant adjustment

self-contained clean water system, single bottle 1.5 liter

quick disconnect for water bottle

wet/dry foot control

Height adjustable work surface, 5"

Central vacuum accessory package asepsis 3-way air/water syringe w/tubing, 5'

autoclavable standard HVE lever valve w/tubing, 5' autoclavable saliva ejector lever valve w/tubing, 5'

Assistant's quad holder assembly Solid work surface, 19" x 24"

Universal heavy duty mounting plate

11" x 14  $^{1}/_{2}$ ", beneath counter (7000-999-BC) Wall mount panel 16" x 18" (7000-999-WM) Air pressure regulator with filter Master air shut-off valve

Umbilical, 7'

# Designer Friendly® Surface Treatments

Powder coat paint color:

Standard colors: O Forest Grey O White

**OR** upgrade powder coat paint (add 10%):

○ Arctic Silver ○ Bronze

O Charcoal

O Copper

Midnight

O Pewter

O Shadow

Solid work surface colors:

**Corian®** 

Antarctica

Gibraltar® Meganite®

○ Arctic Melange
○ Baja Melange

Cinder Granite

O Goldrush Gemstone

O Platinum Granite O Rocky Road Granite

○ Sanibel Granite ○ Moonstone

○ Silver Mist

O Starry Starry Night

Select One Of The Following Three Mounting Options For The 7000INT (N/C):

O Beneath Counter Mount :7000-999-BC

O Under Cabinet Mount :7000-999-UC

O Wall Mount: 7000-999-WM

Cı	ustomize Delivery:	Model Selected:	Part No	Price	Total
Step 1:	Select model and mounting position	Mounting positions for:	7000INT	\$8,270 _	
Step 2:	Select surface treatment options		7000INT-DURR	\$8,770 _	
Step 3: Step 4:	Select handpiece positions Select optional upgrades	Beneath counter mount:	7000-999-BC	NC _	
отер 4.	delect optional apgrades	Under cabinet mount:	7000-999-UC	NC _	
		NC Wall mount:	7000-999-WM	NC	

## Handpiece Positions:

Handpiece Positions	Description	1	2	3	4	5	Part No	Price	Total
SYRINGE	silicone sheath, grey	~					INCLUDED	N/C _	
4-HOLE (vinyl)	vinyl tubing, grey						0017-010-GRY	\$49 _	
■ 6-PIN (silicone)	silicone tubing, grey						0017-614-INT	\$250 _	
SCALER	Satelec SP Newtron Piezo scaler						SAT-S-INT	\$2,615 _	
LED SCALER	Satelec LED Piezo scaler						SAT-S-LED-INT	\$3,355 _	
s additional HP, wet	REQUIRED for a WET 4th HP position						1100-189-INT	\$540 _	
s ADDITIONAL HP, DRY	REQUIRED for a DRY 4th HP position						1100-190-INT	\$270 _	
NO TUBING	leave the block open						OPEN	_	
PLUG BLOCK	plug & close the block						PLUG	_	

<sup>■</sup> REQUIRES HP illumination circuit board & <u>EITHER</u> 250 watt power supply or 25 watt power supply (board can accommodate (2) HP positions) See Page 40 s When a fourth HP position is ordered, an ADDITIONAL DRY or WET HP position must be ordered

## Optional Upgrades:

Description	Part No	Price	Total
Upgrade Kit: wet/dry foot control with air chip	1106-002-INT-CS	\$100 _	
Quick disconnect w/flow control, 1/4" (for water)	0016-454-CS	\$190 _	
Quick disconnect w/flow control, 3/8" (for air)	0016-456-CS	\$190 _	
HP illumination circuit board with time delay and adjustable voltage, accommodates 2 HP positions*			
*Required with 6-pin tubing selection	0014-100	\$280 _	
Power supplies (page 39)			
25 watt power supply (HP illumination)	0014-160	\$160 _	
48 watt power supply (motors)	0014-399	\$325 _	
250 watt universal power supply	0014-387	\$680 _	
Storage organizer with frame and mounting towers	1110-185-CS	\$795 _	
Swivel tray	1110-909-CS	\$645 _	
Additional standard assistant's quad holder bar	0024-156	\$260 _	
Touchpad for 3900 chair (25 ft. cable)	1700-081-1-CS	\$650 _	
Chrome accents	9000-998-CS	\$180 _	
Monitor mount assembly (monitor screen not included)	1300-025-CS	\$570 _	
Standard utility center housing and cover	1115-020	\$190 _	
Additional umbilical (per foot)	1115-999-CS	\$45 _	
Routing valve kit w/water shutoff & regulator/factory installed	1115-315-CS	\$385 _	
		Total: _	

F	Part No	Price
CHAIR ADAPTERS FOREST, 3900 FIXED Accommodates external umbilical	1115-195	\$565
MONITOR POSTS 2" diameter with oilite bushing (used with 1300TP! 30" w/dual slot, fixed sidebox mount	M) 1331-998	\$275
44" FIXED TWO-PIECE LIGHT POST WITH SLOT	OR CORD) 1110-617 1110-616	\$365 \$365

Part No Price



Factory option SAT-S-INT \$2,615

SATELEC SP NEWTRON LED LIT PIEZO SCALER

Factory option SAT-S-LED-INT \$3,355

Please call 800-289-6367 for warranty and technical support of the Scaler.

- · Color coding system ensures proper tip and power setting
- Three tips: #1, #10Z, H3 (Intro Scaling Kit)
- Tips available for periodontics, endodontics, scaling, restorative, and cavity preparation
- Cruise Control™ automatically adjusts frequency amplitude
- On/off controlled by rheostat
- First ever autoclavable LED lit handpiece 100,000 LUX
- · Newtron technology means power, control and reliability!
- Satelec warranty 1 year / tips 30 days from date of purchase
- · Requires additional purchase of Power Supply



#### HANDPIECE ILLUMINATION TUBING, 6 PIN

6-Pin (silicone) tubing, grey

0017-614-INT \$250

- ISO standard power optic handpiece terminal with stainless steel adapter and nut (ISO-L, 6-pin)
- Asepsis style handpiece tubing, 7' (4' to 5' factory installed)



HP ILLUMINATION BOARD WITTH TIME DELAY & VARIABLE VOLTAGE

Light Board 0014-100 \$280

Light board is equipped with a time delay (0-15 seconds) and variable voltages (2.8 - 3.9 VDC) for each handpiece position Power supply not included



## 25 WATT POWER SUPPLY

Factory option

· Required for HP Illumination

0014-160 \$160



#### 48 WATT POWER SUPPLY

Factory option

0014-399

Part No

\$325

Price

- Required for HP Illumination
- Size: 51/2"x 31/2"x1". Easily fits in all Forest floor boxes
- · Quick disconnect connector for Forest scaler and electric motor



#### 250 WATT (UNIVERSAL POWER SUPPLY)

Factory option

0014-387

\$680

- Required for Bien-Air and when using two or more ancillary devices, or when using a
  Forest Operatory light and any Forest ancillary product
  (excluding the 3900 soft-touch massager power source)
- Provides power for MULTIPLE electrical accessories such as: fiber-optics, operatory lights, scalers and electric micro motors
- Size: 71/2" x 51/2" x 21/2". Easily fits in all Forest floor boxes
- Quick disconnect connectors for all of Forest ancillary products
- Built in air-electric switch, enables power source to turn on/off with pilot air
- · Available in both 230 volt and 110 volt



### **ANCILLARY SHELF**

Factory option

0010-310

\$175

- · Accommodates additional ancillary products
- · Designed for easy access to installed devices
- · Powder coated aluminum
- · Installed on IC control head

Part No

Price



HIGH DENSITY WATER BOTTLE		
1 1/2 Liter	0010-464	\$44
1 Liter	0010-364	\$32
1 1/2 Liter with quick disconnect adapter	1100-431	\$93
1 1/2 Liter with quick disconnect upgrade kit	1100-423	\$140



Durr amalgam seperator for fixed chair with sidebox. Cover and mounting bracket included. Factory option 1400-270-CS \$500 Seperator not included. Mounting rack set up for Durr CSI/CASI sperator.



MANUAL SHUTOFF 1100-080 \$38

- Outlet, 3/8" inside diameter compression fitting
- Inlet 1/2" FTP
- · Chrome finished



CONTOURED UTILITY CENTER HOUSING AND COVER (8" x 17 3/16" x 17 3/16"				
White	3900-101	\$555		
Forest Grey	3900-101-G	\$555		
Shadow	3900-101-SH	\$555		
Designer Friendly®	3900-101-C	\$675		
COMPLETE CONTOURED UTILITY CENTER				
<ul> <li>Air and water pressure regulator with filter</li> </ul>				
<ul> <li>Master air and water shut-off valves</li> </ul>				
White	3900-102	\$1,005		
Forest Grey	3900-102-G	\$1,005		
Shadow	3900-102-SH	\$1,005		
Designer Friendly®	3900-102-C	\$1,130		

	Part No	Price
DURR UNIVERSAL HVE VALVE Factory option	1100-505	\$75
DURR SALIVA EJECTOR (SE) VALVE Factory option	1100-510	\$75
DURR SE/ HVE PACKAGE KIT WITH DURR TUBING AND VACUUM Factory option	CANNISTER ( 1100-525	CONNECTORS \$380
DURR HVE VACUUM TUBING (18 MM) Factory option	0017-840	\$90/ meter
DURR HVE VACUUM TUBING (18 MM) Factory option	0017-841	\$85/ meter

		<b>.</b> .	
Part	No	Price	ļ



DURR VACUUM HANGER		
For assemblies:	1091INT, 1092INT	
Factory option	0024-205	\$775
Black	0024-205-BLK	\$775
For assemblies:	1578INT	
Factory option	0024-191	\$775
Black	0024-191-BLK	\$775
For assemblies:	1591INT	
Factory option	0024-215	\$775
Black	0024-215-BLK	\$775



STANDARD DELUXE AIR/WATER SYRINGE WITH QUICK-RELEASE TIP 1109-010 \$115 with sheathed tubing assembly, 5' Grey Silicone 1109-786-INT \$235 STRAIGHT DELUXE AIR/WATER SYRINGE WITH QUICK-RELEASE TIP 1109-025 \$145 with sheathed tubing assembly, 4' Grey Silicone 1109-595-INT \$275



DURR ELECTRIC SWITCH
Factory option 0010-634 \$65





SWIVEL TRAY HOLDER ASSEMBLY (USE WITH ALL STANDARD IC DELIVERIES)
Small swivel tray
1110-468
440
Large swivel tray
1110-475
\$480

Small tray: 213mm x 318mmLarge tray: 289mm x 375mm


Dimensions:


# Pivot Delivery Systems






# Fixed with Sidebox Delivery Systems







# **Utility Centers**




# 3900 Chair




# Chair Mount Light and Monitor Mounts




### 7000-999-BC Beneath Counter




1300-025-CS Monitor Assembly Option




## 7000-999-BC Mounting Plates




0010-086 Heavy Duty Plate




0010-084-CS Option




### 7000-999UC Under Cabinet Mount




### 7000-999-WM Wall Mount




## 7000-999-UC Under Cabinet Leveling/Mounting Plate



## 7000-999-WM Mounting Plate



## 7000-999-WM Leveling/Mounting Board



## Ordering:

Deviation from our standard colors and Designer Friendly® options are subject to an upcharge.

When ordering Forest products, the following information should be included:

**Quantity Accessories Or Options** 

Model/Part #
Date Required
Description

Purchase Order Number

Color

Method Of Shipping

### Work Surface Colors:

Standard Solid Work Surface Colors

(Special non-stock solid work surface colors: call customer service for pricing and quantity discounts.)

Corian® Antarctica

Gibraltar® Arctic Melange, Baja Melange

Meganite® Cinder Granite, Goldrush Gemstone, Moonstone, Platinum Granite, Rocky Road Granite, Sanibel Granite, Silver Mist, Starry Starry Night

## Standard Upholstery Colors (Pearlized colors tend to wear similarly to soft leather)

\*Ultraleather® Apricot, Baltic, Cave, Charcoal, Diplomat Blue, Doe, Dove Grey, Eucalyptus, Lemon, Papyrus, Parrot, Raven Wing, Red, Riviera, Stone. Pearlized Colors: Abalone, Midnight, Pewter \*Ultraleather® is a trademark of Ultrafabrics® LLC.

## Designer Friendly® Paint Colors

Arctic Silver, Bronze, Charcoal, Copper, Grey, Midnight, Pewter, Shadow, White

## Designer Friendly® Plastic Colors

Solid Colors: Forest Grey, White

Carbon Fiber Patterns: Greywood, Lemans, Laser Leather

Tubing Colors: Grey

## Chrome Accents (\$135 upcharge)

Custom Colors deviating from standard subject to a 10%-15% upcharge and extended delivery time

Model Number /Page		Part Number /Page	
Model #	Page	Part #	Page
1091INT-DURR	21	0010-086	48
1092INT-DURR	23	0010-310	40
1095INT	21, 44	0010-364	41
1300RM-INT	25, 27, 28, 47	0010-464	41
1300TPM-INT	25, 47	0014-100	39
1576INT	23, 45	0014-160	40
1577INT	23, 45	0014-381	40
1578INT-DURR	23	0014-387	40
1590INT	21, 44	0014-389	40
1591INT-DURR	21	0017-614-INT	39
2001INT	14, 16, 17, 27, 29, 31, 32, 44	0024-191	43
2003INT	15, 16, 17, 28, 30, 33, 44	0024-205	43
3900INT	13, 46	0024-215	43
3901INT	13, 27, 28, 29, 30, 31, 32, 33, 46	1100-423	41
7000INT	35, 36, 37, 48, 49	1100-505	42
7000INT-DURR	36	1100-510	42
8500INT	22, 27, 28, 29, 30, 45	1100-525	42
8601INT	20, 31, 32, 33, 44	1106-109-INT-CS	41
9073INT	24, 47	1109-010	43
9078INT	24, 28, 32, 47	1109-025	43
9083INT	24, 27, 29, 30, 31, 33, 47	1109-595-INT	43
3907-999-DC	13, 27	1109-786-INT	43
3907-MC-WB	13, 27	1110-468	29, 31, 32, 43
7000-999-BC	35, 36, 37, 48	1110-475	43
7000-999-UC	35, 36, 37, 49	1110-774	30, 33
7000-999-WM	35, 36, 37, 49	1110-787	28
	33, 33, 31, 13	1115-020	46
		1400-270-CS	41
		1700-175-INT	39
		3900-101	41, 46
		SAT-S-INT	39
		SAT-S-LED-INT	39
		OAT O LED INT	00


# 12 Classic




# 5 Pearlized




# Carbon Fiber Colors





# Work Surface Colors



## Warranty

DENTALEZ and its employees are proud of the products we provide to the dental community. We stand behind these products with a warranty against defects in material and workmanship as provided below.

In the event that you experience difficulty with the application or operation of any of our products, please contact our technical service department at our expense at (800)423-3555.

If we cannot resolve the issue by telephone, we will arrange for a representative to contact you or suggest that the product be repaired by a dealer service technician.

If product return is required, we will provide you with a Return Authorization number and shipping instructions to return the product to the proper facility. If the product is under warranty we will ask you to provide proof of purchase such as a copy of your invoice. Please be sure to include the Return Authorization number on the package you are returning. It is the responsibility of the purchaser/end user to immediately cease using the product upon discovery that the product is not in proper working order and to dispose of any defective product that is not returned to DENTALEZ.

Freight costs for product returns are the responsibility of the customer. Products under warranty will be repaired or replaced, at our sole discretion, and returned at our expense. We are not responsible for shipping damages. We will, however, help you file a claim with the freight carrier.

DENTALEZ warrants all equipment and parts to be free of defects in material and workmanship, under normal usage, under the following terms:

#### **Forest Limited Warranty Period\*** Patient chairs / Provider stools 5 years from date of installation Chair / Stool upholstery 2 years from date of installation 5 years from date of installation Delivery systems 7000 series delivery system chassis 20 years from date of installation Handpiece illumination systems 6 months from date of installation Power suppliers, Syringe buttons, HVE and SE valves 1 year from date of installation Water bottles, gas springs, and handpiece tubing 2 years from date of installation Air and Water tubing connectors (T's), Distribution block, 20 years from date of installation Master control block, Flow valves Operatory lights 5 years from date of installation

#### Please note the following additional terms of our warranty and return policy:

- This warranty covers manufacturing defects only and does not cover defects resulting from abuse, accidents, misuse, improper handling, cleaning, care or maintenance. Damage, including, without limitation, discoloration of upholstery, resulting from contact, such as dye transfer, from garments worn by dental staff/team, the use of chemicals, cleaners, disinfectants or sterilization is not covered under this warranty nor is normal wear and tear. Failure to follow or observe our operating/use, maintenance and/or installation instructions voids this warranty as does the use of non-authorized parts and repairs made by an unauthorized repair facility.
- Liability is limited to repair or replacement of the defective part(s) or products at our sole discretion. All other liabilities, in particular liability for damages, including, without limitation, special, exemplary, consequential or incidental damages are excluded.
- THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. NO EMPLOYEE, REPRESENTATIVE OR DEALER IS AUTHORIZED TO CHANGE THIS WARRANTY IN ANY WAY OR TO GRANT ANY OTHER WARRANTY.

#### **CANCELLATION POLICY:**

Standard orders: Written cancellation request required within 3 days of order entry date, otherwise a 25% restock fee will be charged. Expedited ship date requests are non-cancellable.

#### **CUSTOM ORDERS:**

- Custom configurations are those that deviate from Forest's published product offering; these orders are non-cancellable and non-returnable.
- May be subject to an additional "Deviation" charge and/or a non-refundable engineering charge.
- Custom powder coat, custom plastic surface treatments and/or custom treatments such as multiple color combinations or custom colors may result in an upcharge of 10-15% per model and are non-cancellable and non-returnable.

#### **WARRANTY REPLACEMENTS:**

Parts replaced on a product that is in warranty will be warranted for the duration of that product's original warranty

#### PRODUCT RETURN:

Opened products or product returns more than a year old cannot be returned for credit. There will be a 25% (\$25.00 minimum) restocking charge on all items authorized for return.

### NOTES:

\*See individual Product Manuals for complete warranty details. Provided conditions defined in instruction manual are met.

Ultraleather is a registered trademark of Ultrafabrics, Inc., Formica is a registered trademark of The Diller Corporation, Corian is a registered trademark of DuPont, Wilsonart is a registered trademark of Wilsonart LLC., Meganite is a registered trademark of Lai, James, Naugahyde is a registered trademark of United States Rubber Products, Inc., and Active tilt is a registered trademark of Crown Seating.

Designer Friendly is a registered trademark of DTE Oregon, Inc.

### DENTALEZ CORPORATE

6200 NE Cherry Drive Hillsboro, OR 97124

T: 503.640.3012 F: 503.693.9715 DENTALEZ.COM

For all orders and/or questions please call our Customer Experience Representatives at 866.DTE.INFO.

2022 DENTALEZ, Inc. All rights reserved. DENTALEZ, Star and Dentoform are registered trademarks and Nevin is a trademark of DENTALEZ INC. Ramvac is a registered trademark of Ramvac Dental Products Inc.